

Oral Restorative Prosthetics

Authorized 3M™ ESPE™ Lava™ Milling Center

Date Due:

Date Sent: / /

Laboratory: _____

Street: _____

City, State, Zip _____

Contact: _____

Phone #: _____

Dr. Zip: _____ (3M ESPE Required)

Patient: _____

Pontic Design Instructions:

Additional Instructions:

Coping Parameters:

Material: **3M ESPE Lava** [] **3M Essential**: []

Anterior Thickness (.4-2.0 mm): _____

Posterior Thickness (.5-2.0mm): _____

Cement Gap: Tight []
 Default []
 Loose []

Tooth # _____ Shade _____

Tooth # _____ Shade _____

Tooth # _____ Shade _____

Bridge #'s _____ Shade _____

Bridge #'s _____ Shade _____

Finishing of Margins: 95% []
 50% []
 None []
 Other [] _____

Please **do not** mark margins, use any type of die spacer or hardener/sealant on models.

If models are mounted on an articulator, they must be removable from the articulator.

Undercuts, surface defects, occlusal holes, retention grooves should be blocked out with resin block-out.

Severely prepped teeth with missing proximal, buccal or lingual walls must be built-up to adequate prep design

If extensive "wax knife" work is requested or needed for coping build up, additional charges may apply.

Additional model and/or prep work may also incur an additional charge above coping fee.

Please send the following items with this form for the fabrication of Lava Crowns and Bridges.

Pinned and separated models. All sections of the working model must be removable.

For bridges we will need the opposing model and a bite.

Ship case to Oral Restorative Prosthetics, 5200 Snyder Lane, Suite 4, Rohnert Park, CA, 94928

Any questions please give us a call at (707) 584-0775 or fax (707) 584-0705 or email lava@orplab.com

Oral Restorative Prosthetics, 5200 Snyder Lane, Ste. 4, Rohnert Park, CA, 94928